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Anti-citrullinated protein antibodies (ACPA) detection, is the most sensitive and specific test for the diagnosis of Rheumatoid Arthritis (RA). These autoantibodies recognize mainly citrullinated filaggrin, keratin, fibrin, vimentin and collagen and they detected in laboratories using mixtures of cyclic citrullinated peptides (CCP). Anti-CCP antibodies appear several years before the clinical onset of RA and they have higher diagnostic value than rheumatoid factor (RF), antinuclear factor (ANF) and anti-keratin antibodies (AKA). Anti-CCP antibodies are strongly correlated with smoking and common epitope of RA. In addition anti-CCP are positively correlated with disease severity and Larsen score in RA.

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- Suzuki A, Yamada R, Chang X, Tokuhira S, Sawada T, Suzuki M, Nagasaki M, Nakayama-Hamada M, Kawaida R, Ono M, Ohtsuki M, Furukawa H, Yoshino S, Yukioka M, Tohma S, Matsubara T, Wakitani S, Teshima R, Nishioka Y, Sekine A, Iida A, Takahashi A, Tsunoda T, Nakamura Y, Yamamoto K, Functional haplotypes of PADI4, encoding citrullinating enzyme peptidylarginine deiminase 4, are associated with rheumatoid arthritis. *Nat Genet* 2003; 34:395-402.
- Hill JA, Southwood S, Sette A, Levanikar AM, Bell DA, Cairns E, Cutting edge: the conversion of arginine to citrulline allows for a high-affinity peptide interaction with the rheumatoid arthritis-associated HLA-DRB1*0401 MHC class II molecule. *J Immunol* 2003; 171:538-41.
- Lard LR, van Gaalen FA, Schonkeren JJ, Pieterman EJ, Stoeken G, Vos K, Nelissen RG, Westendorp RG, Hoeben RC, Breedveld FC, Toes RE, Huizinga TW, Association of the -2849 interleukin-10 promoter polymorphism with autoantibody production and joint destruction in rheumatoid arthritis. *Arthritis Rheum* 2003; 48:1841-8.
- van Gaalen FA, Visser H, Huizinga TW, A comparison of the diagnostic accuracy and prognostic value of the first and second anti-cyclic citrullinated peptides (CCP1 and CCP2) autoantibody tests for rheumatoid arthritis. *Ann Rheum Dis* 2005; 64:1510-2.
- Wu R, Shovman O, Zhang V, Gilburd B, Zandman-Goddard G, Shoenfeld Y, Increased prevalence of anti-third generation cyclic citrullinated peptide antibodies in patients with rheumatoid arthritis and CREST syndrome. *Clin Rev Allergy Immunol* 2007; 32:47-56.
- Wiik AS, van Venrooij WJ, Pruijn GJ, All you wanted to know about anti-CCP but were afraid to ask. *Autoimmun Rev*.
- Greiner A, Plischke H, Kellner H, Gruber R, Association of anti-cyclic citrullinated peptide antibodies, anti-citrullin antibodies, and IgM and IgA rheumatoid factors with serological parameters of disease activity in rheumatoid arthritis. *Ann N Y Acad Sci* 2005; 1050:295-303.
- Zhao J, Liu X, Wang Z, Liu R, Li Z. Is It Necessary to Combine Detection of Anticitrullinated Protein Antibodies in the Diagnosis of Rheumatoid Arthritis? *J Rheumatol* 2010; in press.
- Rantapaa-Dahlqvist S, de Jong BA, Berglin E, Hallmans G, Wadell G, Stenlund H, Sundin U, van Venrooij WJ, Antibodies against cyclic citrullinated peptide and IgA rheumatoid factor predict the development of rheumatoid arthritis. *Arthritis Rheum* 2003; 48:2741-9.
- Markatseli TE, Voulgari PV, Alamanos Y, Drosos AA, Prognostic Factors of Radiological Damage in Rheumatoid Arthritis: A 10-year Retrospective Study. *J Rheumatol* 2010; in press.
- del Val del Amo N, Ibanez Bosch R, Fito Manteca C, Gutierrez Polo R, Loza Cortina E, Anti-cyclic citrullinated peptide antibody in rheumatoid arthritis: relation with disease aggressiveness. *Clin Exp Rheumatol* 2006; 24:281-6.
- Atzeni F, Turiel M, Capsoni F, Doria A, Meroni P, Sarzi-Puttini P, Autoimmunity and anti-TNF-alpha agents. *Ann N Y Acad Sci* 2005; 1051:559-69.
- Atzeni F, Sarzi-Puttini P, Dell'Acqua D, de Poitu S, Cecchini G, Cruini C, Carrabba M, Meroni PL, Adalimumab clinical efficacy is associated with rheumatoid factor and anti-cyclic citrullinated peptide antibody titer reduction: a one-year prospective study. *Arthritis Res Ther* 2006; 8:R3.
- Braun-Moscovici Y, Markovits D, Zinder O, Schapira D, Rozin A, Ehrenburg M, Dain L, Hoffer E, Nahir AM, Balbir-Gurman A, Anti-cyclic citrullinated protein antibodies as a predictor of response to anti-tumor necrosis factor-alpha therapy in patients with rheumatoid arthritis. *J Rheumatol* 2006; 33:497-500.
- Caro-Oleas JL, Fernandez-Suarez A, Reneses Cesteros S, Porrino C, Nunez-Roldan A, Wichmann Schlipf I, Evaluation of third generation anti-CCP antibodies in the diagnosis of rheumatoid arthritis from undifferentiated polyarthritis after 4 years of follow-up. *Clin Exp Rheumatol* 2008; 26:461-3.
- Quinn MA, Cough AK, Green MJ, Devlin J, Hensor EM, Greenstein A, Fraser A, Emeiy P, Anti-CCP antibodies measured at disease onset help identify seronegative rheumatoid arthritis and predict radiological and functional outcome. *Rheumatology (Oxford)* 2006; 45:478-80.
- Huizinga TW, Amos CI, van der Helm-van Mil AH, Chen W, van Gaalen FA, Jawaheer D, Schreuder GM, Wener M, Breedveld FC, Ahmad N, Lum RF, de Vries RR, Gregersen PK, Toes RE, Criswell LA, Refining the complex rheumatoid arthritis phenotype based on specificity of the HLA-DRB1 shared epitope for antibodies to citrullinated proteins. *Arthritis Rheum* 2005; 52:3433-8.
- Berglin E, Padyukov L, Sundin U, Hallmans G, Stenlund H, Van Venrooij WJ, Klareskog L, Dahlqvist SR, A combination of autoantibodies to cyclic citrullinated peptide (CCP) and HLA-DRB1 locus antigens is strongly associated with future onset of rheumatoid arthritis. *Arthritis Res Ther* 2004; 6:R303-8.
- van Gaalen FA, van Aken J, Huizinga TW, Schreuder GM, Breedveld FC, Zanelli E, van Venrooij WJ, Verweij CL, Toes RE, de Vries RR, Association between HLA class II genes and autoantibodies to cyclic citrullinated peptides (CCPs) influences the severity of rheumatoid arthritis. *Arthritis Rheum* 2004; 50:2113-21.
- van der Helm-van Mil AH, Verpooit KN, Breedveld FC, Huizinga TW, Toes RE, de Vries RR, The HLA-DRB1 shared epitope alleles are primarily a risk factor for anti-cyclic citrullinated peptide antibodies and are not an independent risk factor for development of rheumatoid arthritis. *Arthritis Rheum* 2006; 54:1117-21.
- Matthews AG, Li J, He C, Ott J, Andrade M, Adjusting for HLA-DRbeta1 in a genome-wide association analysis of rheumatoid arthritis and related biomarkers. *BMC Proc* 2009; 3 Suppl 7:S12.
- Klareskog L, Stolt P, Lundberg K, Kallberg H, Bengtsson C, Grunewald J, Ronnelid I, Harris HE, Ulfgren AK, Rantapaa-Dahlqvist S, Eklund A, Padyukov L, Alfredsson L, A new model for an etiology of rheumatoid arthritis: smoking may trigger HLA-DR (shared epitope)-restricted immune reactions to autoantigens modified by citrullination. *Arthritis Rheum* 2006; 54:38-46.